

PLANNING YOUR SPACE

Will it *fit*?

Everything you need to measure a room for the Atlas One and its two support units.

+

≈ 8 m²

+

FLOOR AREA TO ALLOW

+

1.8 m

+

MIN CEILING HEIGHT

+

3

+

UNITS IN TOTAL

+

700 kg

+

CHAMBER WEIGHT

THE THREE UNITS

One chamber plus two support machines that sit nearby and connect to it. Sizes below are the actual units; footprint is the floor space each one takes up.

887 mm

2000

plan

The chamber

Where you sit.

FOOTPRINT 887 × 2000 mm

IN 0.89 × 2.00 m

METRES

HEIGHT 1600 mm (1.60 m)

WEIGHT ≈ 700 kg

660 mm

420

plan

Oxygen & pressure unit

Drives the chamber.

FOOTPRINT 660 × 420 mm

IN 0.66 × 0.42 m

METRES

HEIGHT 900 mm (0.90 m)

WEIGHT ≈ 69 kg

590 mm

390

plan

Air cooler

Keeps it cool.

FOOTPRINT 590 × 390 mm


IN 0.59 × 0.39 m

METRES

HEIGHT 605 mm (0.60 m)


WEIGHT ≈ 43.5 kg

An example layout seen from above, drawn to scale. The chamber sits against a wall with the machines in a row beside it. Allow clear space at the door to get in and out. This is one way to arrange it; the units can be grouped or split to suit your room.



PLAN VIEW · DRAWN TO SCALE · EXAMPLE ARRANGEMENT

The same units shown in three dimensions, so you can picture the heights. The chamber is the large unit; the others are much smaller.



GOOD TO KNOW BEFORE YOU FIT OUT

CEILING	Allow at least 1.8 m of ceiling height. The chamber stands 1.60 m tall and needs clear space above it.	GETTING IN	The chamber weighs about 700 kg and arrives in a large crate. Check doorway widths, lift or stair access, and floor strength on the delivery route.
DOOR SPACE	Keep about 1 m clear at the chamber door so you can get in and out easily.	AIRFLOW	Leave a small gap behind each machine (around 200 mm) for airflow and the connecting hoses.
POWER	Several power points are needed near the machines. Have a licensed electrician confirm dedicated circuits.	FLOOR	A firm, level floor is best. The units are on castors and can be rolled into position.

Please note: the floor plan and clearances above are an example to help you plan. Final placement depends on your room shape, door position, and where power and ventilation are. Confirm the layout on site before fit-out.